

## Product update: **HEINE® LoupeLight 2**

With HEINE HR or HRP Binocular Loupes for a perfect view.

The combination of light and loupe delivers the ultimate in performance as high quality illumination with LED<sup>HQ</sup> maximises the image quality of the loupe.

For a perfect view and a precise diagnosis, the LoupeLight2 is now 25% brighter than its predecessor model, with a homogeneous edge-to-edge illumination of typ. 55,000 lux.

For intensive use, the LoupeLight2 now features a replaceable connecting cord which connects the device to the mPack mini – the new and compact power source for the LoupeLight2 (product description see overleaf).

**NEW!**

- ✓ Now 25% brighter due to improved LED technology
- ✓ Replaceable connecting cord



**HR 2.5x High Resolution Binocular Loupes.**  
Achromatic optics for crisp, bright images.

**Replaceable connecting cord** with high-quality splashproof threaded connector for heavy duty use.

**Now even 25% brighter**

**LED<sup>HQ</sup>**

LED NOW IN HEINE QUALITY.

## Product update: **HEINE® MicroLight 2**

Ultralight, universal headlight.

For those who need a compact, lightweight illumination solution without magnification, we offer the MicroLight 2. Featuring LED<sup>HQ</sup> illumination for homogeneous illumination with true and natural colour rendering for every examination situation. Two mounting options offer individuality and easy handling. Illumination without compromises, in a lightweight design.



HEINE MicroLight 2  
on S-FRAME



HEINE MicroLight 2  
on Lightweight Headband

See overleaf for ordering information.



# HEINE® mPack mini – The compact and light power source for LoupeLight 2 and MicroLight 2.

The new, high performance power source HEINE mPack mini for LoupeLight 2 and MicroLight 2 offers freedom from cables that tie the users to desk, wall or static power sources. With a weight of just 95 g and dimensions of 44 x 23 x 101 mm, this mobile power source is now approximately approx. 70 % smaller and lighter than its predecessor.

The pouch which is included in the delivery can be attached to the belt and the mPack mini can also be charged while clipped to the belt. Thanks to Li-ion battery technology the power source has an operating time of typ. 9 hours, in continuous operation at maximum power typ. 4 hours. The mPack mini can be charged flexibly with USB plug, e.g. on a laptop/desktop or via USB charger.



 USB charging option

**NEW!**  
 ✓ 70 % smaller and lighter  
 ✓ Operating time typ. 9 hours  
 ✓ Charging time typ. 4 hours



Size comparison of HEINE mPack LL/mPack mini

## LoupeLight2 Sets with HR Binocular Loupes

HR 2.5x optics with i-View on S-FRAME, LoupeLight 2, detachable yellow filter, mPack mini, mPack mini Belt Clip, E4-USB plug-in transformer, protective lenses, 2 sterilisable levers, retaining cord, cleaning fluid, cleaning cloth, carrying case

340 mm working distance	C-008.32.450
420 mm working distance	C-008.32.451
520 mm working distance	C-008.32.452

## MicroLight 2 Sets

MicroLight 2 on S-FRAME with mPack mini, mPack mini Belt Clip, E4-USB plug-in transformer, retaining cord and cleaning fluid in case

MicroLight 2 on Lightweight Headband with mPack mini, mPack mini Belt Clip, E4-USB plug-in transformer, case

## Single Devices

LoupeLight 2 for HR/HRP Loupes with i-View loupe mount	C-008.32.226
mPack mini with Li-ion rechargeable battery and E4-USB plug-in transformer	X-007.99.650
mPack mini with Li-ion rechargeable battery, without E4-USB plug-in transformer	X-007.99.649
Connecting cord SC 1 (1.6 m) for LoupeLight 2 or MicroLight 2 / mPack mini	C-000.32.240



Имя HEINE характеризует продукцию высочайшего качества!

**HEINE** QUALITY  
MADE IN GERMANY

"

-

"

<http://heine-med.ru/>  
[info@heine-med.ru](mailto:info@heine-med.ru)

: (495) 902-59-26 . / ; (495) 518-55-99.

HEINE,

HEINE ( ) -

!

