



∴ HEINE BETA[®] 4 Rechargeable Handles HEINE[®] NT 4 Table Charger

Next-generation energy management



∴ HEINE BETA® 4 Rechargeable Handles HEINE® NT 4 Table Charger

Product Summary

Charge your HEINE BETA 4 USB handle with Li-ion technology with any USB plug or other USB power source. The BETA 4 USB handle has an integrated, exclusive Galvanic Separation Module that protects the handle, and the user, even if a faulty or inferior plug is used; thus preventing the possibility of the full mains voltage flowing directly to the handle. A charge status indicator always lets you know the charge status of your handle. Compatible with all HEINE 3.5 V instrument heads.



PERFORMANCE
FUNCTION
DURABILITY



∴ HEINE BETA® 4 Rechargeable Handles

Performance Function Durability



Feature

4 handle options:

- ∴ BETA dry cell handle
- ∴ BETA 4 USB rechargeable handle with medical approved E 4 USB plug-in power supply
- ∴ BETA NT 4 and BETA 4 SLIM NT rechargeable handle for desk charging in the HEINE NT 4 table charger



Benefit

- ∴ Complete freedom of choice
- ∴ Numerous applications



∴ HEINE BETA® 4 Rechargeable Handles

Performance Function Durability



Feature

Li-ion technology

Benefit

- ∴ No “memory effect”
- ∴ Charge your handle how you like without worrying about charge level or battery condition
- ∴ Higher capacity



∴ HEINE BETA® 4 Rechargeable Handles

Performance Function Durability



Feature

Deep Battery Discharge Protection

Benefit

- ∴ Prevents batteries from completely draining
- ∴ Ensures maximum battery life
- ∴ Protects your investment



∴ HEINE BETA® 4 Rechargeable Handles

Performance Function Durability



Feature

BETA 4 Handles are compatible with all 3.5 V XHL & LED HEINE instruments

Benefit

Compatibility during your entire career



∴ HEINE BETA® 4 Rechargeable Handles

Performance Function Durability



- ∴ Slightly knurled handle finish: non-slip grip, easy to clean and disinfect
- ∴ High quality metal construction: strength, durability
- ∴ Exclusive Automatic Connector with locking design: no metal on metal friction, secure fit for a lifetime
- ∴ 5 years guarantee





USB



∴ HEINE BETA[®] 4 USB Handle



∴ HEINE BETA® 4 USB Rechargeable Handle

Performance Function Durability



Feature

USB charging

Benefit

New flexible and convenient charging option:

- ∴ Laptop or computer
- ∴ Medical approved E4 USB plug-in power supply
- ∴ Other USB chargers




USB



∴ HEINE BETA® 4 USB Rechargeable Handle

Performance Function Durability



Feature

Integrated, exclusive Galvanic Separation Module

Benefit

- ∴ Prevents the possibility of the full mains voltage to flow directly to the handle
- ∴ Protects the handle and the user, even if a faulty or inferior plug is used



USB



∴ HEINE BETA® 4 USB Rechargeable Handle

Performance Function Durability



Feature

HEINE BETA4 USB Bottom Insert

- ∴ Charge Status Indicator



USB

Benefit

- ∴ During charging process: indicator blinks green
- ∴ Rechargeable battery is fully charged: indicator shines green permanently
- ∴ Battery must be recharged: indicator shines orange



HEINE QUALITY
MADE IN GERMANY



∴ HEINE® NT 4 Table Charger



∴ HEINE® NT 4 Table Charger

Performance Function Durability



Feature

Quick Charge

Benefit

Approx. 2-6 hours to recharge
(with BETA 4 NT or BETA 4 SLIM NT
rechargeable handle)



∴ HEINE® NT 4 Table Charger

Performance Function Durability

Feature

HEINE NT4 Table Charger:

- ∴ Exclusive Ready-to-use Function
- ∴ Handles switch ON & OFF automatically when placed into charger
- ∴ Once your desired brightness level is set, it stays there every time you use your instrument



Benefit

Prevents battery damage and premature bulb burn outs, if the instruments are accidentally left on





**∴ HEINE® BETA 4 NT and BETA 4 SLIM NT
Rechargeable Handles**



∴ HEINE BETA® 4 NT Rechargeable Handles

Performance Function Durability

Feature



HEINE BETA 4 NT Bottom Insert:

- ∴ Charge Status Indicator
- ∴ During charging process: the illuminated ring of the NT4 pulses blue
- ∴ Rechargeable battery is fully charged: illuminated ring of the NT4 shines permanently
- ∴ Battery must be recharged: BETA4NT Rechargeable Handle light shines orange

Benefit

Never run out of battery power unexpectedly during examination



∴ HEINE® BETA 4 NT and BETA 4 SLIM NT Rechargeable Handles Essentials



NEW!

- ∴ Li-ion technology
- ∴ Exclusive Automatic Connector
- ∴ Charge Status Indicator
- ∴ USB – the new, flexible charging option

