


Instrument Head for HEINE UniSpec® Tubes



[01]



[02]



[03]

Illumination by Fiber Optic cable and projector, lamp handle or 3.5V BETA Handle with illumination adaptor (X-002.98.002).

- ∴ **Annular Fiber Optic illumination.** Bright, reflex-free illumination.
- ∴ **Magnifying lens 1.5x.** High-resolution image.
- ∴ **Hinged viewing window with insufflation port.** Improved vision by insufflation.
- ∴ **Stainless Steel.**

Instrument head for UniSpec Disposable Tubes	
complete with viewing window and swivel lens 1,5x [01]	E-003.18.100
without viewing window and swivel lens	E-003.18.098
Viewing window only [02]	E-003.18.099
Swivel lens 1.5x, only (for sigmoidoscopes) [03]	E-003.18.096
same, but 1.8x (only suitable for proctoscopes and anoscopes)	E-003.19.099

HEINE UniSpec® Disposable Tubes



85x20 mm 130x20 mm 250x20 mm

- ∴ **Disposable.** No cross-infection.
- ∴ **No cleaning required.** Saves time.
- ∴ **Moulded in patient-friendly thermoplastic.**
No cold sensation on insertion.
- ∴ **Light transmission through the tube wall.** Distal annular light exit.
Bright, homogeneous illumination of the area examined, reflex reduction.
- ∴ **Graduated inside the tube for patient comfort.** Indication of insertion length.

UniSpec Disposable Anoscopes (85x20 mm)	
Pack of 25 disposable anoscopes	E-003.19.925
Box of 100 (4 packs of 25) anoscopes	E-003.19.911

UniSpec Disposable Proctoscopes (130x20 mm)	
Pack of 25 disposable proctoscopes	E-003.19.825
Box of 100 (4 packs of 25) proctoscopes	E-003.19.811

UniSpec Disposable Sigmoidoscopes (250x20 mm)	
Pack of 25 disposable sigmoidoscopes	E-003.18.825
Box of 100 (4 packs of 25) sigmoidoscopes	E-003.18.811

⊗: Single use only.



HEINE Anoscope/Proctoscope Illumination Head



[01]

One single illumination head for all the tubes shown on this page.

For 2.5V or 3.5V handle or Fiber Optic illumination with Fiber Optic projector and cable.

- ∴ **A wide range of tubes with either direct or 90° direction of view available.** Many different diagnostic and treatment methods are possible.
- ∴ **Tubes can be rotated through 360°.** Select the ideal position for treatment or examination.
- ∴ **Viewing window and swivel lens both offer 2x magnification.** High-resolution image.
- ∴ **Smooth, rounded surface.** Easy to clean. F.O. version can be cleaned with conventional disinfectants.



[02]



[03]

Illumination head	2.5V XHL	3.5V XHL	F.O.
complete with swivel lens and viewing window [01]	E-001.19.101	E-002.19.101	E-003.19.120
Illumination head without accessories	E-001.19.102	E-002.19.102	E-003.19.121
XHL Xenon Halogen spare bulb	X-001.88.040	X-002.88.051	
Swivel lens, only [02]			E-000.19.103
Viewing window, only [03]			E-000.19.104

Tubes with straight direction of view

Anoscope/Proctoscope with proximal illumination.

- ∴ **Distal end of tube angled.** Assures a good view and easy access with instruments.
- ∴ **All-metal construction.** Easy to clean and disinfect with conventional solutions.



60x22mm

120x22mm

Tubes with straight direction of view

Anoscope tube 60mm x 22mm dia., complete with obturator	E-000.19.111
Proctoscope tube 120mm x 22mm dia., complete with obturator	E-000.19.107

Tubes with 90° direction of view (Blond)

Proctoscope with proximal illumination.

- ∴ **Polished metal viewing mirror.** Tough, unbreakable.
- ∴ **End-piece with mirror detachable.** Easy to clean, disinfect or sterilize.



100x20mm



130x25mm

Proctoscope Tubes

Proctoscope tube, 90° direction of view, 130mm x 25mm dia.	E-000.19.105
Proctoscope tube, 90° direction of view, 100mm x 20mm dia.	E-000.19.114

HEINE F.O. Anoscopes with annular Fiber Optic Illumination



70x25mm 70x20mm 70x15mm

Illumination by Fiber Optic cable and projector, lamp handle or 3.5V BETA handle with illumination adaptor (X-002.98.002).

- ∴ **Three different diameters: 15 mm, 20 mm and 25 mm.** Useful range of sizes.
- ∴ **Annular Fiber Optic illumination.** Bright illumination of the area to be examined.
- ∴ **Angled tube-end.** Ensures a good view and easy access with instruments.
- ∴ **Rounded proximal shield.** Easy insertion and rotation.
- ∴ **Tube can be dismantled.** Easy to clean and disinfect with conventional solutions.
- ∴ **Stainless steel.** Durable.

Anoscope	F.O.
70x25 mm dia., complete with obturator	E-003.19.325
70x20 mm dia., complete with obturator	E-003.19.320
70x15 mm dia., complete with obturator	E-003.19.315

HEINE F.O. Proctoscopes with Distal Annular Illumination



Illumination by Fiber Optic cable and projector, lamp handle or 3.5V BETA handle with illumination adaptor (X-002.98.002).

- ∴ **Distal annular Fiber Optic (F.O.) illumination.** Bright, reflex-free illumination.
- ∴ **Magnifying lens 1.8x magnification.** High-resolution image.
- ∴ **Hinged viewing window with insufflation port.** Vision can be improved by insufflation.
- ∴ **Angled tube-end.** Improves vision and access.
- ∴ **End-piece can be unscrewed.** Easy to clean and disinfect with conventional solutions.
- ∴ **Stainless steel.** Durable.

Proctoscope	F.O.
130x20 mm dia., complete with obturator, viewing window, swivel lens and end-piece	E-003.19.213
Proctoscope tube without accessories [01]	E-003.19.013
Obturator, only [02]	E-003.19.113
End-piece, only [03]	E-003.19.095
Viewing window, only [04]	E-003.18.099
Swivel lens 1.8x, only [05]	E-003.19.099

HEINE F.O. Sigmoidoscopes with Distal Annular Illumination

Illumination by Fiber Optic cable and projector, lamp handle or 3.5V BETA handle with illumination adaptor (X-002.98.002).

- ∴ Available in various sizes: 300x20 mm, 250x20 mm, 200x20 mm, 250x16 mm, 200x11.4 mm. Sizes suitable for most patients of all ages.
- ∴ Distal annular Fiber Optic (F.O.) illumination. Bright, reflex-free illumination.
- ∴ Swivel lens with 1.5x magnification. High-resolution image.
- ∴ Opening viewing window with insufflation port. Vision can be improved by insufflation.
- ∴ Graduated in cm. Indication of insertion length.
- ∴ Plastic end-piece. Insulated.
- ∴ End-piece can be unscrewed. Easy to clean and disinfect with conventional solution.
- ∴ Stainless steel. Durable.



Sigmoidoscope	F.O.
300x20 mm dia. with obturator, viewing window, swivel lens and end-piece [01]	E-003.18.230
End-piece, only	E-003.18.095
250x20 mm dia. with obturator, viewing window, swivel lens and end-piece	E-003.18.225
End-piece, only	E-003.18.095
200x20 mm dia. with obturator, viewing window, swivel lens and end-piece	E-003.18.220
End-piece, only	E-003.18.095
250x16 mm dia. with obturator, viewing window, swivel lens and end-piece [02]	E-003.18.525
End-piece, only	E-003.18.088
200x11.4 mm dia. with obturator, viewing window and swivel lens [03]	E-003.18.240
End-piece, only	E-003.18.094
Viewing window, only [04]	E-003.18.099
Swivel lens 1.5x, only [05]	E-003.18.096

Telescope



For use with proctoscope and rectoscope with distal illumination and with instrument head for UniSpec disposable tubes.

- ∴ **Superb prismatic optical system with colour-corrected lenses.** Upright image in true colour.
- ∴ **Magnification 3 – 5.5x.** Excellent resolution.
- ∴ **Hinged clip-fitting.** Does not hinder instrumentation.

Telescope

E-000.18.908

Rubber Insufflation Bulb/Hygiene Filter



- ∴ **Rubber insufflation** bulb to use with anoscopies, proctoscopes and sigmoidoscopes.
- ∴ **Hygiene filter** to use with the rubber insufflation bulb. This filter is intended for single use only and prevents contamination of the rubber insufflation bulb.

Rubber insufflation bulb

E-000.18.105

Disposable Hygiene Filters 25 pcs/pack

E-000.18.116

Sponge Holder



- ∴ **High quality design in stainless steel.**

Sponge holder 400 mm working length

E-000.18.906

HEINE Illumination Adaptor



Can be attached to any HEINE 3.5V power source to illuminate any HEINE fiber optics instrument, e.g. anoscope, proctoscope or sigmoidoscope.

Illumination adaptor

3.5V XHL

Illumination adaptor

X-002.98.002

XHL Xenon Halogen spare bulb

X-002.88.078

Anoscope Instrument Kit ATE

An attractively-priced kit, recommended for diagnosis and therapy of haemorrhoids.



Set complete with: **Instrument head** for UniSpec Disposable Tubes with viewing window and swivel lens 1.5x.
Pack of 25 UniSpec Anoscope Tubes 85 mm
Plug-in transformer E7 with lamp handle

Anoscope instrument kit ATE	6V XHL
Complete Anoscope instrument kit ATE	E-095.16.510
XHL Xenon Halogen spare bulb	X-004.88.060
Lamp Handle with switch and XHL Xenon Halogen bulb	X-004.99.607
Cord switch (suitable for HEINE instruments and power sources with cinch connector)	X-095.16.202
Plug-in transformer E7	X-095.16.104

Sigmoidoscope/Proctoscope Kit RTE

Inexpensive but complete kit of instruments for routine rectal examination.



Set complete with: **Instrument head** for UniSpec Disposable Tubes with viewing window and swivel lens 1.5x.
Pack of 25 UniSpec Sigmoidoscope Tubes 250 mm
Pack of 25 UniSpec Proctoscope Tubes 130 mm
Pack of 25 UniSpec Anoscope Tubes 85 mm
Sponge holder
Insufflation bulb
Plug-in transformer E7 with lamp handle

Sigmoidoscope/Proctoscope kit RTE	6V XHL
Complete Sigmoidoscope/Proctoscope kit RTE	E-095.16.501
Disposable Hygiene Filters 25 pcs/pack	E-000.18.116
XHL Xenon Halogen spare bulb	X-004.88.060
Lamp Handle with switch and XHL Xenon Halogen bulb	X-004.99.607
Cord switch (suitable for HEINE instruments and power sources with cinch connector)	X-095.16.202
Plug-in transformer E7	X-095.16.104

Sigmoidoscope/Proctoscope Kit RE 7000

This kit is similar to RTE but features even better illumination by Fiber Optic Projector HK7000.



Set complete with: **Instrument head** for UniSpec Disposable Tubes with viewing window and swivel lens 1.5x

Pack of 25 UniSpec Sigmoidoscope Tubes 250mm

Pack of 25 UniSpec Proctoscope Tubes 130mm

Pack of 25 UniSpec Anoscope Tubes 85mm

Sponge holder

Insufflation bulb

F.O. Projector HK7000, 150W

Fiber Optic cable

Sigmoidoscope/Proctoscope Kit RE 7000	F.O.
Complete Sigmoidoscope/Proctoscope kit	E-096.15.501
Spare bulb, 150W	Y-096.15.102
Disposable Hygiene Filters 25 pcs/pack	E-000.18.116

Имя HEINE характеризует продукцию высочайшего качества!

HEINE QUALITY
MADE IN GERMANY

"

-

"

<http://heine-med.ru/>
info@heine-med.ru

: (495) 902-59-26 .. , (495) 518-55-99

HEINE,

HEINE () -

!

